

File :Q:\GCMS-08 Minimal\DATA\20221213\05- -jop.D  
Operator : Organic QC Lab  
Acquired : 13 Dec 2022 12:43 using AcqMethod Hy\_VOC-SV\_24cm\_S100.M  
Instrument : GCMS 8  
Sample Name: EPH-WA-ALI-R1  
Misc Info : aliphatics in dcm  
Vial Number: 4



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 2:1, 2µl inj.; GC/MS; Injector 270°C

