

File :L:\QC Data\GCMS-04 Minimal\DATA\010720\06-.D  
 Operator :  
 Acquired : 7 Jan 2020 13:46 using AcqMethod VOC-200SPLIT DECKER.M  
 Instrument : GCMS4  
 Sample Name: FIA-CAL-07  
 Misc Info :  
 Vial Number: 70



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

