

File : Q:\GCMS-04 Minimal\DATA\080823A\35 GA.D
Operator :
Acquired : 9 Aug 2023 16:07 using AcqMethod VOC-100 SPLIT #H2 28MIN(HOT).M
Instrument : GCMS4
Sample Name: FIA-OLE
Misc Info :
Vial Number: 16



Column: DB-624UI, 20m x 0.18mm ID x 1.0µm
Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold
GC Parameters: 1µl inj.; GC-MS; Injector temp 250°C carrier gas Hydrogen

