

File :Q:\GCMS-08 Minimal\DATA\20220228B\22- =JP=.D
Operator : Organic QC Lab
Acquired : 28 Feb 2022 19:48 using AcqMethod Hy_BDE_24cm_S1.M
Instrument : GCMS 8
Sample Name: FRS-023S
Misc Info : Tris(2,3-Dibromopropyl)phosphate; in toluene
Vial Number: 59



Column: DB-5, 30m, 0.25 ID, 0.25um
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj; GC/MS; INJ 250c.

