

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-04 Minimal\DATA\1219823\
... 37- KM.D
Operator :
Instrument : GCMS4
Acquired : 20 Dec 2023 9:38 using AcqMethod VOC-20 SPLIT#H2 28MIN255DEG HOT.M
Sample Name: FTOH-028S
Misc Info : 1H,1H,13H-Perfluorotridecan-1-ol in MeOH



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

