

File : M:\GCMS-07 Minimal\DATA\08112101\49-ALR.D
Operator : Organic QC department
Acquired : 12 Aug 2021 01:35 using AcqMethod KERLA-2019_Short_15Min_2S_2uL.M
Instrument : 5977B MSD
Sample Name: H-001S
Misc Info : 1-Methylnaphthalenein Toluene
Vial Number: 35



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

