

File :Q:\GCMS-06 Minimal\DATA\011921\18- KM.D
Operator : Organic QC Lab
Acquired : 19 Jan 2021 13:37 using AcqMethod CHICK_2017_13MIN_S100.M
Instrument : GCMS 6
Sample Name: HS-003S-D-40X
Misc Info : RT: Naphtha in MeCl2
Vial Number: 35



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

