

File :D:\MassHunter\GCMS\1\data\12121901\11-.D  
Operator : Organic QC department  
Acquired : 12 Dec 19 01:44 pm using AcqMethod KERLA-2Split-2uL.M  
Instrument : 5977B MSD  
Sample Name: ISO6468-BENZ  
Misc Info :in Hexane  
Vial Number: 53



Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
Oven Program: 80c, 17c/min to 340c, 8min  
GC Parameters: Cons. Split, 12psi constant flow  
Split 20:1, 1uL inj.; GC/MS; INJ 270c

