

File : L:\QC Data\GCMS-08 Minimal\DATA\20220614\11-1 - SA
... K.D
Operator : Organic QC Lab
Instrument : GCMS 8
Acquired : 14 Jun 2022 15:46 using AcqMethod Hy_Pest_15cm_S20.M
Sample Name: L-010N
Misc Info :



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

