

File :Q:\GCMS-04 Minimal\DATA\032323A\37-RM.D
Operator :
Acquired : 24 Mar 2023 7:22 using AcqMethod VOC-100 SPLIT #H2 20MIN175DEG.M
Instrument : GCMS4
Sample Name: M-502-47-10X
Misc Info :
Vial Number: 59



Column: DB-624UI, 20m x 0.18mm ID x 1.0µm
Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold
GC Parameters: 1µl inj.; GC-MS; Injector temp 250°C carrier gas Hydrogen

