

File : Q:\GCMS-04 Minimal\DATA\1219823\21 RM.D  
 Operator :  
 Acquired : 20 Dec 2023 00:44 using AcqMethod VOC-100 SPLIT #H2 20MIN175DEG.M  
 Instrument : GCMS4  
 Sample Name: M-502A-R  
 Misc Info :  
 Vial Number: 21



Column: DB-624UI, 20m x 0.18mm ID x 1.0µm  
 Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold  
 GC Parameters: 1µl inj.; GC-MS; Injector temp 250°C carrier gas Hydrogen

