

File :L:\QC Data\GCMS-07 Minimal\DATA\10132202\31 -~~SA
... K.D
Operator : Organic QC department
Instrument : 5977B MSD
Acquired : 14 Oct 2022 00:44 using AcqMethod KERLA-2Split-2uL_VOC-SV_W-200EM.M
Sample Name: M-552.2-LPC
Misc Info :



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5m in, 10psi constant flow
Split 2:1, 2µl inj.; GC/MS; Injector 270°C

