

File : Q:\GCMS-07\DATA\03291901\02- D  
 Operator : K.D.  
 Acquired : 29 Mar 19 07:21 am using AcqMethod KERLA-2017.M  
 Instrument : 5977B MSD  
 Sample Name: M-610-QC-FL  
 Misc Info : PAH Performance check in ACCN  
 Vial Number: 16

Column: DB-5, 30m, 0.25 ID, 0.25um  
 Oven Program: 100c, 1m, 30c/min to 320c, final 2min.  
 GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow  
 Split 20:1, 1ul inj.; GC/MS; INJ 250c.

