

File :D:\MassHunter\GCMS\1\data\10291901\05- .D  
 Operator : KD & DC  
 Acquired : 29 Oct 19 09:03 am using AcqMethod KERLA-2Split-2uL.M  
 Instrument : 5977B MSD  
 Sample Name: M-610-QC  
 Misc Info :  
 Vial Number: 83

Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
 Oven Program: 80c, 17c/min to 340c, 8min  
 GC Parameters: Cons. Split, 12psi constant flow  
 Split 20:1, 1uL inj.; GC/MS; INJ 270c

