

File :Q:\GCMS-07 Minimal\DATA\05042101\41- -DC.D
Operator : Organic QC department
Acquired : 04 May 2021 23:27 using AcqMethod KERLA-2017_Short_15Min.M
Instrument : 5977B MSD
Sample Name:M-625-18
Misc Info : Phenol-D5 in DCM
Vial Number: 93



Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

