

File :M:\GCMS-06 Minimal\DATA\081321\02- ALR.D
 Operator : Organic QC Lab
 Acquired : 13 Aug 2021 14:29 using AcqMethod CHICK_2019_VOC-SV_S100.M
 Instrument : GCMS 6
 Sample Name: M-625-BN
 Misc Info : Base/Neutral Extractables in DCM
 Vial Number: 25



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
 Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
 GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
 Split 20:1, 1µl inj.; GC/MS; Injector 270°C

