

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-04 Minimal\DATA\062322\3  
 ... RM.D  
 Operator :  
 Instrument : GCMS4  
 Acquired : 23 Jun 2022 14:43 using AcqMethod VOC-100SPLIT#.M  
 Sample Name: M-8010B  
 Misc Info :



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

