

File :K:\QC Data\GCMS-04 Minimal\DATA\082120A\07-220081297-RM.D
Operator : ORGANIC QC DEPT.
Acquired : 21 Aug 2020 17:36 using AcqMethod VOC-100SPLIT HOT_235 D.M
Instrument : GCMS4
Sample Name: M-8015B-5031-24 (220081297)
Misc Info :
Vial Number: 22



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

