

File :Q:\GCMS-04 Minimal\DATA\052423\50- (JP).D  
Operator :  
Acquired : 25 May 23 4:06 pm using AcqMethod VOC-20 SPLIT #H2 20MIN175DEG.M  
Instrument : GCMS4  
Sample Name: M-8032 (223051599)  
Misc Info : Acrylamide; in meoh  
Vial Number: 41



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

