

File : Q:\GCMS-07\DATA\09201801\11-
 Operator : K.D
 Acquired : 20 Sep 2018 13:08 using AcqMethod KERLA-2017.M
 Instrument : 5977B MSD
 Sample Name: M-8060
 Misc Info: in Isooctane
 Vial Number: 7

Column: DB-5, 30m, 0.25 ID, 0.25um
 Oven Program: 100c, 1m, 30c/min to 320c, final 2min.
 GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow
 Split 20:1, 1ul inj.; GC/MS; INJ 250c.

