

File :D:\MassHunter\GCMS\1\data\04291901\16-
Operator : KD
Acquired : 29 Apr 2019 10:17 using AcqMethod KERLA-2017.M
Instrument : 5977B MSD
Sample Name: M-8141A-IC
Misc Info : in Hexane
Vial Number: 43

Column: DB-5, 30m, 0.25 ID, 0.25um
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

