

File : Q:\GCMS-08 Minimal\DATA\20231030\11- (JP).D
Operator : Organic QC Lab
Acquired : 30 Oct 2023 18:21 using AcqMethod Hy_Pest_15cm_S5.M
Instrument : GCMS 8
Sample Name: M-8150-02-PFB
Misc Info : 2,4-D-PFB; in mtbe
Vial Number: 6



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

