

File :D:\MassHunter\GCMS\1\data\01141901\04-  
 Operator : K.D  
 Acquired : 14 Jan 2019 09:09 using AcqMethod KERLA-2017.M  
 Instrument : 5977B MSD  
 Sample Name: M-8270-10-ASL  
 Misc Info : in CH2C12  
 Vial Number: 3

Column: DB-5, 30m, 0.25 ID, 0.25um  
 Oven Program: 100c, 1m, 30c/min to 320c, final 2min.  
 GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow  
 Split 20:1, 1ul inj.; GC/MS; INJ 250c.

