

File : L:\QC Data\GCMS-06 Minimal\DATA\062923\65 -SAK.D
 Operator : Organic QC Lab
 Acquired : 30 Jun 2023 06:05 using AcqMethod CHICK_2019.M
 Instrument : GCMS 6
 Sample Name: M-GAC
 Misc Info :
 Vial Number: 31

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
 Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
 GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
 Split 20:1, 1µl inj.; GC/MS; Injector 270°C

