

File :Q:\GCMS-04\DATA\021419\45-.D
 Operator : RFXM
 Acquired : 15 Feb 19 10:36 am using AcqMethod VOC-200SPLIT DECKER.M
 Instrument : GCMS4
 Sample Name: M-GRA-CAL-R-IS-R-05
 Misc Info :
 Vial Number: 39



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

