

File :Q:\GCMS-04 Minimal\DATA\101821B\59 YX.D
Operator : Organic QC Lab
Acquired : 20 Oct 2021 10:08 using AcqMethod VOC-200SPLIT HOT260.M
Instrument : GCMS4
Sample Name: M-GRA-VT-IS-AS
Misc Info : D5580 Valve Timing Cablibration
Vial Number: 12



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

