

File :Q:\GCMS-04 Minimal\DATA\032823B\79 RM.D
 Operator :
 Acquired : 29 Mar 2023 23:02 using AcqMethod VOC24-100SPLIT#H2 25MIN 255DEG.M
 Instrument : GCMS4
 Sample Name: M-GRO-EPA-CC
 Misc Info :
 Vial Number: 44



Column: DB-624UI, 20m x 0.18mm ID x 1.0µm
 Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold
 GC Parameters: 1µl inj.; GC-MS; Injector temp 250°C carrier gas Hydrogen

