

File : Q:\GCMS-04 Minimal\DATA\110521\15- RM.D

Operator :

Acquired : 5 Nov 2021 15:07 using AcqMethod VOC-SCAN 24 10MIN.M

Instrument : GCMS4

Sample Name: M-PAPCLUS

Misc Info :

Vial Number: 31

(221111114)

Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm

Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)

GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

