

File :Q:\GCMS-04 Minimal\DATA\122122\17- RM.D  
Operator :  
Acquired : 21 Dec 2022 20:25 using AcqMethod VOC-100 SPLIT #H2 5MIN.M  
Instrument : GCMS4  
Sample Name: M-REF-X-09  
Misc Info :  
Vial Number: 37



Column: DB-624UI, 20m x 0.18mm ID x 1.0µm  
Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold  
GC Parameters: 1µl inj.; GC-MS; Injector temp 250°C carrier gas Hydrogen

