

File :Q:\GCMS-06 Minimal\DATA\090221\BSB\24- KM.D
Operator : Organic QC Lab
Acquired : 02 Sep 2021 16:55 using AcqMethod CHICK_2019_S100.M
Instrument : GCMS 6
Sample Name: MO-COMP-D-40X (221091050)
Misc Info : Motor Oil Composite Standard @ 5000 ug/mL in MeCl2
Vial Number: 40



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

