

File :L:\QC Data\GCMS-07 Minimal\DATA\02162201\30 -SAK.D
Operator : Organic QC department
Acquired : 16 Feb 2022 18:17 using AcqMethod KERLA-2017_Short_15Min.M
Instrument : 5977B MSD
Sample Name: P-002S
Misc Info : Aldrin in MeOH
Vial Number: 149



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

