

**AccuStandard®**

**Sample Name: P-017R-WL-10X in iso**

File: C:\CPData\DATA\20240425\224041471 JP(1).0000.RAW

Date (start): 4/25/2024 5:32:58 PM

Column: 110c, 0min, 15c/min to 325, Final 2 min

Parameters: Splitless; 2ul; ECD: 325c; Inj: 250c

notes:

Method: Ni\_FAST\_ECD\_S10.MET


Setpoint File: Ni\_FAST\_ECD\_S10.M68

Calibration File:

Inst: GC #7

Oper: JJP

Vial = 2



10.94	Tetrachloro-meta-xylene (877-09-8)
16.37	Chlordane (Technical Mixture) (12789-03-6)
22.57	Decachlorobiphenyl (2051-24-3)