

Operator : Organic OC Lab  
File : D:\MassHunter\GCMS\1\data\20200729\57- JP.D  
Acquired : 30 Jul 2020 12:46  
Instrument : GCMS 8  
Sample Name: P-170S using AcqMethod Hy\_Pest\_15cm\_S20.M



Column: DB-5, 30m, 0.25ID, 0.25um  
Oven: 80c, 2m, 30c/m to 320, final 2min  
GC para: Pulse split 40psi 0.5m, 10psi constant flow  
Split 20:1, 1ul inj, GCMS, Inj 280c

Misc Info : Phorate in meoh  
Vial Number: 66

