

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-09 Minimal\DATA\02152023
... \36 KM.D
Operator :
Instrument : GCMS 9
Acquired : 15 Feb 2023 20:28 using AcqMethod Hy SV DB-5MS UI20M 0.18 40 split Time 14 MIN .M
Sample Name: P-196S
Misc Info : Dacthal in MeOH



Column: DB-5MS UI, 20m x 0.18mm ID x 0.18µm
Oven Program: 45°C, 13°C/min to 290°C (hold 5min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector temp 270°C

