

File :Q:\GCMS-08 Minimal\DATA\20200123\16- .D  
Operator :  
Acquired : 23 Jan 2020 16:38 using AcqMethod Hy\_Pest\_15cm\_S20.M  
Instrument : GCMS 8  
Sample Name: P-220S-H-10X  
Misc Info : EPN in Hexane  
Vial Number: 31



**AccuStandard®**

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

