

File :D:\MassHunter\GCMS\1\DATA\04252202\22-ALR.D  
Operator : Organic QC department  
Acquired : 25 Apr 2022 18:11 using AcqMethod KERLA-2017.M  
Instrument : 5977B MSD  
Sample Name: P-415S-CN  
Misc Info : Methyl nonyl ketone in AccCN  
Vial Number: 112



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

