

File :Q:\GCMS-06 Minimal\DATA\062320\22-YX.D  
Operator : Organic QC Lab  
Acquired : 23 Jun 2020 22:32 using AcqMethod CHICK\_2019\_13MIN.M  
Instrument : GCMS 6  
Sample Name: P-693S  
Misc Info : Pyridaben in MeOH  
Vial Number: 67



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25 $\mu$ m  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1 $\mu$ l inj.; GC/MS; Injector 270°C

