

File :Q:\GCMS-08 Minimal\DATA\20210707\85-YX.D
Operator : Organic QC Lab
Acquired : 08 Jul 2021 19:06 using AcqMethod Hy_PAH_15cm_S5.M
Instrument : GCMS 8

Misc Info : Trimethyl phosphate (TMP) in toluene
Vial Number: 132



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

