

File :Q:\GCMS-06\DATA\011119\23-.D  
 Operator : D.C  
 Acquired : 12 Jan 19 01:53 am using AcqMethod CHICK\_2017\_S100.M  
 Instrument : GCMS 6  
 Sample Name: PHENOL-QME  
 Misc Info : QME Mix in MeOH  
 Vial Number: 51



Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
 Oven Program: 80c, 17c/min to 340c, 8min  
 GC Parameters: Cons. Split, 12psi constant flow  
 Split 100:1, 1uL inj.; GC/MS; INJ 270c

