

File :Q:\GCMS-08 Minimal\DATA\20211015b\03-JP.D
Operator : Organic QC Lab
Acquired : 15 Oct 2021 14:11 using AcqMethod Hy_PAH_24cm_S50.M
Instrument : GCMS 8
Sample Name: PLAS-PL-103S
Misc Info : t-Butylphenyl diphenyl phosphate; hex; Tech Mix of isomers
Vial Number: 52



Column: DB-5, 30m, 0.25ID, 0.25um
Oven: 80c, 2m, 30c/m to 320, final 2min
GC para: Pulse split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj, GCMS, Inj 280c

