

File :Q:\GCMS-08 Minimal\DATA\20191118\27- .D  
Operator : KM / JP / ZA  
Acquired : 18 Nov 2019 22:28 using AcqMethod KERLA-2017-100uL draw speed.M  
Instrument : GCMS 8  
Sample Name: PLAS-UV-007S-CN  
Misc Info : 2-(5-tert-Butyl-2-hydroxyphenyl)benzotriazole in AcCN  
Vial Number: 32



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

