

File :B:\QC Data\GCMS-08 Minimal\DATA\20190912\46-.D
Operator : JJP / KM /
Acquired : 13 Sep 2019 10:59 using AcqMethod KERLA-2017.M
Instrument : GCMS 8
Sample Name: RAC-05
Misc Info : 4-CHLOROANILINE IN ACCN
Vial Number: 146



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

