

File : Q:\GCMS-08 Minimal\DATA\20200331\15- .D  
Operator : Organic QC Lab  
Acquired : 31 Mar 2020 18:47 using AcqMethod Hy\_PAH\_24cm\_S100.M  
Instrument : GCMS 8  
Sample Name : SK-200X-100ML  
Misc Info : Sulfur in Light distillate kerosene  
Vial Number : 79



Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
Oven Program: 80c, 24c/min to 340c, 2min  
GC Parameters: Cons. Split, 5.4psi constant flow  
Split 100:1 1uL inj.; GC/MS; INJ 270c

