

File :Q:\GCMS-08 Minimal\DATA\20231017\24- (JP).D
 Operator : Organic QC Lab
 Acquired : 17 Oct 2023 17:23 using AcqMethod Hy_PAH_24cm_S200.M
 Instrument : GCMS 8
 Sample Name: Z-014C-R
 Misc Info :in hex:tol [1:1]
 Vial Number: 4



Column: DB-5, 30m, 0.25 ID, 0.25um
 Oven Program: 80c, 1m, 15c/min to 340c, final 2min
 GC Parameters: Split 10:1, 1ul inj.; GC/MS; INJ 270c
 (Actual analytical conditions may vary from these general parameters)

