

File : L:\QC Data\GCMS-06 Minimal\DATA\051322\20 - SAK.D  
Operator : Organic QC Lab  
Acquired : 13 May 2022 17:57 using AcqMethod CHICK\_2017\_13MIN\_S100.M  
Instrument : GCMS 6  
Sample Name : Z-014F-2  
Misc Info :  
Vial Number : 122



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

